



+ NEW PICO[◇] 14 sNPWT

designed to challenge
hard-to-heal wounds

PICO[◇] 14

Single Use Negative Pressure
Wound Therapy System

Helping you get **CLOSER TO ZERO[◇]**
delay in wound healing
smith-nephew.com/pico

Our new **PICO 14 Single Use Negative Pressure Wound Therapy System (sNPWT)** with **AIRLOCK[◇] Technology**

- Pump duration of up to 14 days¹
- Aimed for use on deep wounds greater than 2cm deep²
- Enhanced pump to aid use in large wounds with less user intervention^{3†}

Turn around wound healing trajectory more effectively than standard dressings⁴ and tNPWT⁵ with PICO sNPWT.

Smith+Nephew
















For detailed product information, including indications for use, contraindications, precautions and warnings, please consult the product's applicable Instructions for Use (IFU) prior to use.

[†] Compared to previous versions. **References:** **1.** Smith+Nephew December 2018. PICO 14 Service Life Testing: 14 Day Device Lifespan. Internal Report. RD/18/132. **2.** Smith+Nephew 2018. The review of evidence supporting the use of PICO in wounds ≥ 2 cm in depth. Internal report EO AWM. PCS230.001 v2. **3.** Smith+Nephew 2019. Research & Development Report. PICO 14 and PICO 7 Initial Pump Down and Maintenance Pump Down Time Outs RD/19/084. **4.** Dowsett C, Hampton J, Myers D, Styche T. Use of PICO to improve clinical and economic outcomes in hard-to-heal wounds. *Wounds International*. 2017;8, p53–58. **5.** Kirsner R, Dove C, Reyzelman A, Vayser D, Jaimes H. A prospective, randomized, controlled clinical trial on the efficacy of a single-use negative pressure wound therapy system, compared to traditional negative pressure wound therapy in the treatment of chronic ulcers of the lower extremities. *Wound Repair Regen*. 2019 May 14. [Epub ahead of print].

•Trademark of Smith+Nephew ©August 2019 Smith+Nephew. 21013 | CA98757 03/20

Product ordering codes



	PICO 14 device + 2 dressings	Multipack with 5 dressings	
Dressing sizes	Code	Code	
 Multisite small 15cm x 20cm	66802040	66802020	
 Multisite large 20cm x 25cm	66802041	66802021	
 10cm x 20cm	66802042	66802022	
 15cm x 15cm	66802045	66802025	
 15cm x 20cm	66802046	66802026	
 15cm x 30cm	66802047	66802027	
 20cm x 20cm	66802048	66802028	
 25cm x 25cm	66802049	66802029	
Consumables	Code	Antimicrobial	Code
 Foam dressing filler 10cm x 12.5cm	66801021	 ACTICOAT®Flex 3 10 cm x10 cm	66800399
 RENASYS®gel patch 10cm x 7cm	66801082	 ACTICOAT®Flex 7 2.5cm x 60cm	66800545

*Dressing kit comprising of 2 dressings, 1 pump and pump clip, secondary fixation strips. Low exuding wounds – up to 14 day system. Additional dressings are required for moderately exuding wounds (available separately).

**Multipack of dressings comprising of 5x individually packaged sterile dressings, secondary fixation strips.